



Robert F. Faulkner, D.D.S., Inc.

Maxillofacial Prosthodontist

Biteplane Care Instructions

- ❖ Never bite the biteplane into place. Put it in by hand. If it feels off or is not going into place remove it completely from the mouth and start over.
- ❖ Brush your biteplane twice daily or when you brush your teeth. You can use Softsoap antibacterial hand soap or toothpaste. When not in use, keep hydrated in a container with water. If you are traveling, wrap the biteplane in a moist paper towel and place in your case. Place the case in a ziplock baggie to keep from drying out.
- ❖ Keep out of extreme temperatures to prevent warping.
- ❖ **KEEP OUT OF REACH OF PETS!!** Your pet will chew this like a toy, especially since it already has your scent. Keep on a high counter or in your medicine cabinet.
- ❖ It is very important to bring your biteplane with you to **EVERY** appointment you may have with our office.